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人proBNP、BNP和NT-proBNP

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- Benjamin, EJ et al. Heart disease and stroke statistics –2017 update. *Circulation*. 2017, 135 (10): e146-e603.
- Emdin, M et al. Comparison of brain natriuretic peptide (BNP) and amino-terminal ProBNP for early diagnosis of heart failure. *Clin Chem*. 2007, 53: 1289-1297.
- Mueller, T et al. Diagnostic accuracy of B type natriuretic peptide and amino terminal proBNP in the emergency diagnosis of heart failure. *Heart*. 2005, 91: 606-612.
- Ponikowski, P et al. 2016 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure. *Eur. Heart J*. 2016, 37: 2129-2200.
- Schellenberger, U et al. The precursor to B-type natriuretic peptide is an O-linked glycoprotein. *Arch Biochem Biophys*. 2006, 451: 160-166.
- Semenov, AG et al. Processing of proBNP is suppressed by O-glycosylation in the region close to the cleavage site. *Clin Chem*. 2009, 55(3): 489-498.
- Mair, J et al. The impact of cardiac natriuretic peptide determination on the diagnosis and management of heart failure. *Clin Chem Lab Med*. 2001, 39: 571-588.
- Cowie, MR and Mendez, GF. BNP and congestive heart failure. *Prog Cardiovasc Dis*. 2002, 44: 293-321.
- Pandey, KN. Biology of natriuretic peptides and their receptors. *Peptides*. 2005, 26: 901-932.
- Belenky, A et al. The effect of class-specific protease inhibitors on the stabilization of B-type natriuretic peptide in human plasma. *Clin Chim Acta*. 2004, 340: 163-172.
- Murdoch, DR et al. Disparity between studies of the stability of BNP in blood: comparison of endogenous and exogenous peptide. *Heart*. 1999, 81: 212-213.
- Niederkofer, EE et al. Detection of endogenous B-type natriuretic peptide at very low concentrations in patients with heart failure. *Circ. Heart Fail*. 2008, 1(4): 258-264.

13. Brandt, I et al. Dipeptidyl-peptidase IV converts intact B-type natriuretic peptide into its des-SerPro form. *Clin Chem.* 2006, 52(1): 82-87.
14. Pankow, K et al. Successive action of meprin A and neprilysin catabolizes B-type natriuretic peptide. *Circ Res.* 2007, 101(9): 875-882.
15. Müller, D et al. Rat insulin-degrading enzyme: cleavage pattern of the natriuretic peptide hormones ANP, BNP, and CNP revealed by HPLC and mass spectrometry. *Biochemistry.* 1992, 31(45): 11138-11143.
16. Ralat, LA et al. Insulin-degrading enzyme modulates the natriuretic peptide-mediated signaling response. *J Biol Chem.* 2011, 286(6): 4670-4679.
17. Seferian, KR et al. The brain natriuretic peptide (BNP) precursor is the major immunoreactive form of BNP in patients with heart failure. *Clin Chem.* 2007, 53: 866-873.
18. Luckenbill, KN et al. Cross-reactivity of BNP, NT-proBNP, and proBNP in commercial BNP and NT-proBNP assays: preliminary observations from the IFCC Committee for Standardization of Markers of Cardiac Damage. *Clin Chem.* 2008, 54(3): 619-621.
19. Semenov, AG et al. Searching for a BNP standard: glycosylated proBNP as a common calibrator enables improved comparability of commercial BNP immunoassays. *Clin Biochem.* 2017, 50: 181-185.
20. Saenger, AK et al. Specificity of B-type natriuretic peptide assays: cross-reactivity with different BNP, NT-proBNP, and proBNP peptides. *Clin Chem.* 2017, 63(1): 351-358.
21. Seferian, KR et al. Immunodetection of glycosylated NTproBNP circulating in human blood. *Clin Chem.* 2008, 54: 866-873.
22. Nishikimi, T et al. The effect of glycosylation on plasma N-terminal proBNP-76 levels in patients with heart or renal failure. *Heart.* 2012, 98(2): 152-161.
23. Røsjø, H et al. Influence of glycosylation on diagnostic and prognostic accuracy of N-terminal pro-B-type natriuretic peptide in acute dyspnea: data from the Akershus Cardiac Examination 2 Study. *Clin Chem.* 2015, 61(8): 1087-1097.
24. Giuliani, I et al. Assay for measurement of intact B-type natriuretic peptide prohormone in blood. *Clin Chem.* 2006, 52: 1054-1061.
25. Liang, F et al. Evidence for functional heterogeneity of circulating B-type natriuretic peptide. *J Am Coll Cardiol.* 2007, 49: 1071-1078.
26. Waldo, SW et al. Pro-B-type natriuretic peptide levels in acute decompensated heart failure. *J Am Coll Cardiol.* 2008, 51: 1874-1882.
27. Dries, DL et al. Simultaneous assessment of unprocessed proBNP1-108 in addition to processed BNP32 improves identification of high-risk ambulatory patients with heart failure. *Circ Heart Fail.* 2010, 3(2): 220-227

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Semenov AG, Postnikov AB, Tamm NN, Seferian KR, Karpova NS, Bloshchitsyna MN, Koshkina EV, Krasnoselsky MI, Serebryanaya DV, Katrukha AG. Processing of pro-brain natriuretic peptide is suppressed by O-glycosylation in the region close to the cleavage site. *Clin Chem.* 2009 Mar;55(3):489-598.

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Semenov AG, Katrukha AG. Different Susceptibility of B-Type Natriuretic Peptide (BNP) and BNP Precursor (proBNP) to Cleavage by Neprilysin: The N-Terminal Part Does Matter. Clin Chem. 2016 Apr;62(4):617-22

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Semenov AG, Katrukha AG. Analytical Issues with Natriuretic Peptides – has this been Overly Simplified? EJIFCC. 2016 Jul; 27(3): 189–207.

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Semenov AG, Tamm NN, Apple FS, Schulz KM, Love SA, Ler R, Feygina EE, Katrukha AG. Searching for a BNP standard: Glycosylated proBNP as a common calibrator enables improved comparability of commercial BNP immunoassays. Clin Biochem. 2017 Mar;50(4-5):181-185.

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Immunoassay Kit for Quantification of BNP and proBNP (US 9,145,459)

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