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I. Postnikov et al., 2012: Prediction of short-term to medium-term cardiac events and death in patients with suspected acute coronary syndrome

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II. Hjortebjerg et al., 2015: Prediction of long-term cardiac death in patients with Type 1 Diabetes

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		Cardiovascular Mortality	
		HR (95% CI)*	P Value
NT-IGFBP-4	Unadjusted	6.58 (3.75-11.57)	<.001
	Adjusted for nephropathy	4.13 (2.16-7.88)	<.001
	Adjusted for nephropathy and cardiovascular risk factors	5.81 (2.62-12.86)	<.001
CT-IGFBP-4	Unadjusted	4.90 (2.70-8.82)	<.001
	Adjusted for nephropathy	2.82 (1.43-5.55)	.003
	Adjusted for nephropathy and cardiovascular risk factors	2.58 (1.10-6.10)	.030

*HR: hazard ratio; CI: confidence interval

III. Hjortebjerg et al., 2017: Prediction of short-term to long-term mortality in patients with acute myocardial infarction

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IV. Postnikov et al., 2018: Prediction of short-term to medium-term all-cause mortality in patients with acute heart failure

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1. Postnikov, A. B. et al. N-terminal and C-terminal fragments of IGFBP-4 as novel biomarkers for short-term risk assessment of major adverse cardiac events in patients presenting with ischemia. Clin. Biochem. 45, 519-524 (2012).
2. Hjortebjerg, R. et al. PAPP-A and IGFBP-4 fragment levels in patients with ST-elevation myocardial infarction treated with heparin and PCI. Clin. Biochem. 48, 322-328 (2015).
3. Hjortebjerg, R. et al. IGFBP-4 Fragments as Markers of Cardiovascular Mortality in Type 1 Diabetes Patients With and Without Nephropathy. The Journal of Clinical Endocrinology & Metabolism 100, 3032-3040 (2015).
4. Konev, A. A. et al. Characterization of endogenously circulating IGFBP-4 fragments—Novel biomarkers for cardiac risk assessment. Clinical Biochemistry 48, 774-780 (2015).
5. Konev, A. A. et al. Glycosylated and non-glycosylated NT-IGFBP-4 in circulation of acute coronary syndrome patients. Clin. Biochem. (2018). doi:10.1016/j.clinbiochem.2018.03.004
6. Postnikov, A. B. et al. C-terminal fragment of IGFBP-4 is independently associated with mortality in patients hospitalized due to acute heart failure. Poster presentation at ESC Heart Failure 2018, <http://u.to/MmpyEg> (P1828).
7. Hjortebjerg, R. et al. Insulin-Like Growth Factor Binding Protein 4 Fragments Provide Incremental Prognostic Information on Cardiovascular Events in Patients With ST-Segment Elevation Myocardial Infarction. J Am Heart Assoc 6, (2017).

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Detection of IGFBP-4 Fragments as a Diagnostic Method

EP 2448969, US 9012610, US 9964549, JP 5840605, KR 101767611 and CN 102471379. This invention describes a method for diagnosing of cardiovascular diseases, which comprises detection of IGFBP-4 (insulin-like growth factor binding protein-4) fragments in patients' blood. It provides antibodies as well as epitopes for antibodies, specific to proteolytic fragments (both N- and C-terminal) of IGFBP-4 originated from IGFBP-4 molecule after its cleavage by specific protease PAPP-A. Antibodies could be used for development of immunoassay methods for quantitative or qualitative detection of IGFBP-4 fragments in human blood.

Method for Determining the Risk of cardiovascular Events Using IGFBP Fragments

US 10191066. This invention describes the method for determining the risk of future major adverse cardiovascular events, which comprises detection of proteolytic fragments of IGFBP-4 or IGFBP-5 (insulin-like growth factor binding protein 4 or insulin-like growth factor binding protein 5) in patients' blood. The present invention provides antibodies and immunoassays, suitable for specific measurement of proteolytic fragments of IGFBPs. In current invention the IGFBP fragments are suggested to be utilized as blood biomarkers for the risk prediction of major adverse cardiovascular events (MACE).

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