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1994

2014

Together. Today and Tomorrow.

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Filatov VL, et al. Troponin: structure, properties, and mechanism of functioning. *Biochemistry.* 1999, 64(9): 969-985.

Katrakha A, et al. Biochemical factors influencing measurement of cardiac troponin I in serum. *Clin. Chem. Lab. Med.* 1999, 37(11-12): 1091-1095.

Katrakha A, et al. New approach to standardization of human cardiac troponin I (cTnI). *Scand. J. Clin. Lab. Invest., Suppl.* 1999, 230: 124-127.

Katrakha AG, et al. A new method of human cardiac troponin I and troponin T purification. *Biochem. Mol. Biol. Int.* 1995, 36: 195-202.

Katrakha AG, et al. Troponin I is released in bloodstream of patients with acute myocardial infarction not in free form but as complex. *Clin. Chem.* 1997, 43(8): 1379-1385.

Katrakha AG, et al. Degradation of cardiac troponin I: implication for reliable immunodetection. *Clin. Chem.* 1998, 44(12): 2433-2440.

Katrakha AG. Antibody selection strategies in cardiac troponin assays. *Cardiac Markers*, 2nd edition, Edited by Alan H. B. Wu. 2003, 173-185.

Postnikov AB, et al. N-terminal and C-terminal fragments of IGFBP-4 as novel biomarkers for short-term risk assessment of major adverse cardiac events in patients presenting with ischemia. *Clin. Chem.* 2012, 45: 519-524.

Røsjø H, et al. Diagnostic utility of a single-epitope sandwich B-type natriuretic peptide assay in stable coronary artery disease: data from the Akershus Cardiac Examination (ACE) 1 Study. *Clin. Biochem.* 2012, 45(16-17): 1269-1275.

Schulz O, et al. Clinical differences between total PAPP-A and measurements specific for the products of free PAPP-A activity in patients with stable cardiovascular disease. *Clin. Biochem.*

Seferian KR, et al. The brain natriuretic peptide (BNP) precursor is the major immunoreactive form of BNP in patients with heart failure. *Clin. Chem.* 2007, 53: 866-873.

Seferian KR, et al. Immunodetection of glycosylated NT-proBNP circulating in human blood. *Clin. Chem.* 2008, 54(5): 866-873.

Semenov AG and Seferian KR. Biochemistry of the human B-type natriuretic peptide precursor and molecular aspects of its processing. *Clin. Chim. Acta.* 2011, 412(11-12): 850-860.

Semenov AG, et al. Processing of pro-brain natriuretic peptide is suppressed by O-glycosylation in the region close to the cleavage site. *Clin. Chem.* 2009, 55(3): 489-498.

Semenov AG, et al. Processing of Pro-B-Type Natriuretic Peptide: Furin and Corin as Candidate Convertases. *Clin. Chem.* 2010, 56(7): 1166-1176.

Semenov AG, et al. Human Pro-B-Type Natriuretic Peptide Is Processed in the Circulation in a Rat Model. *Clin. Chem.* 2011, 57(6): 883-890.

Tamm NN, et al. Novel immunoassay for quantification of brain natriuretic peptide and its precursor in human blood. *Clin. Chem.* 2008, 54(9), 1511-1518.

Tamm NN, et al. Measurement of B-type natriuretic peptide by two assays utilizing antibodies with different epitope specificity. *Clin. Biochem.* 2011, 44(2-3): 257-259.

Vylegzhanina AV, et al. Epitope Specificity of Anti-Cardiac Troponin I Monoclonal Antibody 8I-7. *Clin. Chem.* 2013, 59(12): 1814-1816.

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- Alhadi HA and Fox KA.** Do we need additional markers of myocyte necrosis: the potential value of heart fatty-acid-binding protein. *QJMed.* 2004; 97: 187-198.
- Bayes-Genis A. et al.** Pregnancy-associated plasma protein A as a marker of acute coronary syndromes. *N Engl J Med.* 2001, 345(14): 1022-9.
- Body R and Ferguson C.** Pregnancy-associated plasma protein A: a novel cardiac marker with promise. *Emerg Med J.* 2006, 23(11): 875-877.
- Brennan ML et al.** Prognostic value of myeloperoxidase in patients with chest pain. *New Eng J Med.* 2003, 349(17): 1595-1604.
- Clearfield MB.** C-reactive protein: a new risk assessment tool for cardiovascular disease. *JAOA* 2005; 105(9):409-416.
- Hájek P et al.** Pregnancy-associated plasma protein A and proform eosinophilic major basic protein in the detection of different types of coronary artery disease. *Physiol. Res.* 2008, 57(1): 23-32.
- Heeschen C et al.** Pregnancy-associated plasma protein-A levels in patients with acute coronary syndromes Comparison with markers of systemic inflammation, platelet activation, and myocardial necrosis. *J. Am. Coll. Cardiol.* 2005, 45:229-237.
- Katrukha AG.** Antibody selection strategies in cardiac troponin assays. *Cardiac Markers*, 2nd edition, Edited by Alan HB. Wu. 2003, 173-185.
- Katrukha AG et al.** Troponin I is released in bloodstream of patients with acute myocardial infarction not in free form but as complex. *Clin. Chem.* 1997, 43(8): 1379-1385.
- Katrukha AG et al.** Degradation of cardiac troponin I: implication for reliable immunodetection. *Clin. Chem.* 1998, 44(12): 2433-2440.
- Kiss A and Reinhart W.** [Determination of myoglobin in serum and in urine after heart infarction]. *Wien Klin. Wochenschr.* 1956, 68(9):154-155.
- Koenig W et al.** Prospective Study of High-Sensitivity C-Reactive Protein as a Determinant of Mortality: Results from the MONICA/KORA Augsburg Cohort Study, 1984-1998. *Clin. Chem.* 2008; 54:2:335-342.
- Meyer MH et al.** CRP determination based on a novel magnetic biosensor. *Biosens. Bioelectron.* 2007, 22(6):973-979.
- Nambi V.** The use of myeloperoxidase as a risk marker for atherosclerosis. *Curr. Atheroscler. Rep.* 2005, 7(2):127-131.
- Penttilä K et al.** Myoglobin, creatine kinase MB isoforms and creatine kinase MB mass in early diagnosis of myocardial infarction in patients with acute chest pain. *Clin. Biochem.* 2002, 35(8):647-653.
- Postnikov AB et al.** N-terminal and C-terminal fragments of IGFBP-4 as novel biomarkers for short-term risk assessment of major adverse cardiac events in patients presenting with ischemia. *Clin. Chem.* 2012, 45:519-524.
- Qin QP et al.** Release Patterns of Pregnancy Associated Plasma Protein A (PAPP-A) in Patients with Acute Coronary Syndroms. *Scand Cardiovasc. J.* 2002, 36(6):358-361.
- Ridker PM.** C-reactive protein: a simple test to help predict risk of heart attack and stroke. *Circulation.* 2003. 108:e81-e85.
- Scirica BM et al.** Clinical Application of C-Reactive Protein Across the Spectrum of Acute Coronary Syndromes. *Clin. Chem.* 2007, 53:1800-1807.
- Seferian KR et al.** Immunodetection of glycosylated NT-proBNP circulating in human blood. *Clin. Chem.* 2008, 54(5):866-873.
- Shiesh SC et al.** Determination of C-reactive protein with an ultra-sensitivity immunochemiluminometric assay. *J. Immunol. Methods.* 2006, 311(1-2):87-95.
- Sin KK et al.** Fluorogenic nanocrystals for highly sensitive detection of C-reactive protein. *IEE Proc. Nanobiotechnol.* 2006, 153(3):54-58.
- Storch J and Thumser AE.** The fatty acid transport function of fatty acid-binding proteins. *Biochim. Biophys. Acta* 2000; 1486: 28-44
- Terkelsen CJ. et al.** Temporal course of pregnancy-associated plasma protein-A in angioplasty-treated ST-elevation myocardial infarction patients and potential significance of concomitant heparin administration. *Am. J. Cardiol.* 2009, 103(1):29-35.
- Tertti R et al.** Intravenous administration of low molecular weight and unfractionated heparin elicits a rapid increase in serum pregnancy-associated plasma protein A. *Clin Chem.* 2009, 55(6):1214-1217.
- van Nieuwenhoven FA et al.** Discrimination between myocardial and skeletal muscle injury by assessment of the plasma ratio of myoglobin over fatty acid-binding protein. *Circulation.* 1995, 92(10):2848-2854.
- Zhang R. et al.** Association between myeloperoxidase levels and risk of coronary artery disease. *JAMA.* 2001, 286(17):2136-2142.

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