



# DataSheet

**CATALOGUE #:** 4CA72

**PRODUCT NAME:** Monoclonal mouse anti-CA72-4

---

**MAbs:** 1C2, 7C1

Hybridoma clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with native CA72-4.

**MAb isotypes:** IgG1 for MAbs 1C2, 7C1

**Applications:** Detection of CA72-4 antigen.

Recommended pair for sandwich immunoassay (capture – detection): 1C2 – 7C1

MAbs are working in Western blotting.

**Purification:** Ion exchange chromatography.

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)