

CATALOGUE #: RC4TC2

PRODUCT NAME: Recombinant anti-cardiac troponin complex

MAbs:	RecChim20C6 Recombinant full-sized chimeric 20C6 expressed in mammalian cells. Original wild type 20C6 variable domains and human IgG1 constant domains are combined in this construct.
Specificity:	Human cardiac troponin complex
Applications:	Human cardiac troponin complex sandwich immunoassays
Purification:	Protein A chromatography
Presentation:	PBS, pH 7.4, 0.09 % sodium azide (NaN ₃)
Storage:	+4 °C (+2 ... +8 °C allowed)
Material safety note:	<p>This product is sold for research or further manufacturing use only. Standard Laboratory Practices should be followed when handling this material.</p> <p>Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.</p>