



# DataSheet

**CATALOGUE #:** 3C13

**PRODUCT NAME:** Monoclonal mouse anti-cyclosporin A

---

**MAbs:** CSZ22

**Specificity:** Cyclosporin A  
MAb CSZ22 recognizes cyclosporin A in buffer solutions and whole blood extracts.

**MAb isotypes:** **IgG1** for MAb CSZ22

**Purification:** Protein G chromatography

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6  
FI-20520 Turku FINLAND  
[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)