



DataSheet

CATALOGUE #: 3L1

PRODUCT NAME: Monoclonal mouse anti-*Listeria monocytogenes*

MAbs: LZF7, LZH1

Specificity: MAbs recognize *Listeria monocytogenes* outer membrane fraction and intact *L. monocytogenes* cells.

MAb isotypes: **IgG1** for MAb LZH1
IgG2a for MAb LZF7

Applications: Indirect ELISA and Western blotting
Recommended pair for sandwich immunoassay (capture-detection): LZH1 – LZF7.

Purification: Protein G chromatography.

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi