



DataSheet

CATALOGUE #: 4ER1

PRODUCT NAME: Monoclonal mouse anti-erythropoetin

MAbs: Epo1, Epo2

MAb isotypes: IgG1 for MAbs Epo1, Epo2

Purification: Protein G chromatography

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.