



DataSheet

CATALOGUE #: 4IT3

PRODUCT NAME: Monoclonal mouse anti-immunoreactive trypsin (IRT)

MAbs: 2C4, 4A7

Hybridoma clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with highly purified human trypsin.

MAb isotypes: IgG for MAbs 2C4, 4A7

Applications: MAb are working in immunoassays.

Recommended pairs for sandwich immunoassay (capture-detection): 4A7– 2C4

Purification: Ion exchange

Presentation: PBS, 0.05 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note:

This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi