



# DataSheet

**CATALOGUE #:** 5F3cc

**PRODUCT NAME:** Monoclonal mouse anti-FITC

---

**MAbs *in vitro*:** 2A3cc

Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with FITC-BSA conjugate.

**Specificity:** Fluorescein, FITC-labeled proteins.

2A3cc reacts with both free and conjugated FITC.

**MAb isotypes:** IgG1 for 2A3cc

**Applications:** Detection of FITC-labeled antigens or antibodies.

Capture antibody for FITC-labeled biopolymers.

Enzyme-labeled second antibody for immunohistochemistry with FITC-conjugated mouse monoclonal antibodies.

**Purification:** Protein A chromatography

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligent 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)