

# Datasheet

Blood coagulation and Anemia • Bone Metabolism • Cardiac Markers • Fertility and Pregnancy  
Gangliosides • Hormone Markers • Immunology and Serology • Infectious Diseases • Inflammation  
Kidney Diseases • Metabolic Syndrome • Microbial and Plant Toxins • Miscellaneous • Neuroscience  
Thyroid Diseases • Tumor Markers • Veterinary

**CATALOGUE #:** 8CA19

**PRODUCT NAME:** CA19-9

<b>Source:</b>	Human metastatic liver carcinoma. Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.
<b>Applications:</b>	Immunogen, calibrator in immunoassays, in vitro studies.
<b>Specific activity:</b>	Specific activity confirmed by enzyme immunoassay.
<b>Presentation:</b>	PBS, pH 7.4, 3 % sucrose, 0.05 % of sodium azide (NaN <sub>3</sub> ).
<b>Storage:</b>	-20 °C (-15 ... -30 °C allowed)
<b>Material safety note:</b>	This product is sold <b>for research or further manufacturing use only</b> . Standard Laboratory Practices should be followed when handling this material. Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.