



# DataSheet

**CATALOGUE #:** 8G16-2b

**PRODUCT NAME:** Monosialoganglioside GM1, bovine

---

**Source:** Bovine brain tissue.

**Applications:** Useful as an antigen and receptor for cholera toxin, growth inhibitor marker in fibroblasts and marker for lymphoid subpopulations. Stable for the majority of sialidase.

**Analysis:** Purity more than 98 % by HPTLC. MW 1545

**Presentation:** Lyophilized from water solution.

**Storage:** Before and after reconstitution at -20 °C (-15 ... -30 °C allowed)

**Material safety note:** This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6  
FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)