



DataSheet

CATALOGUE #: 8G16-9b

PRODUCT NAME: Disialoganglioside GD3, bovine

Source: Bovine brain tissue.

Applications: One of the most important markers for malignant melanoma. Antibodies to the GD3 ganglioside can induce partial remission of tumor growth in animals as well as in humans via enhancement of cytotoxic and proliferative response of lymphocytes.

Analysis: Purity more than 98 % by HPTLC. MW 1461

Presentation: Lyophilized from water solution.

Storage: Before and after reconstitution at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi