



DataSheet

CATALOGUE #: PCC2

PRODUCT NAME: Polyclonal anti-cystatin C

Immunogen: Recombinant human cystatin C (HyTest Cat.# 8CY5)

Host animal: Sheep

Specificity: Human Cystatin C

Applications: Human Cystatin C immunoassay, immunoprecipitation, immunohistochemistry, Western blotting

Purification: Immunoaffinity method utilizing cystatin C-Sepharose as an affinity matrix.

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi